

Joint Analytical Systems (JAS), based in Moers Germany, is a world leader in providing unique GC solutions into the laboratory and online marketplace. Our leadership position was accomplished by providing economically priced solutions to customers' analytical requirements using Agilent Technologies instruments as a solution platform.

Since May 2002, we have inherited the rights to manufacture, sell, support, and develop the former Agilent Atomic Emission Detector (AED). The AED, a gas chromatography (GC)-based detector, allows detection of virtually all elements, with excellent selectivity, which is widely used in routine QA/QC and research laboratories. **For AED applications we have developed important relationships:**

Diablo Analytical, Inc., an independent company that provides R&D level analytical application development services for the AED and other analytical technologies. Diablo Analytical will continue to provide application support and training for JAS AED customers as they have since 1996 for Agilent Technologies. Diablo Analytical can be contacted at 925-609-1150 or by e-mail info@diabloanalytical.com

weMeasureIt, are experts in analytical techniques, gas chromatographic solutions and software development to enhance product quality, product compliance and process performance. With weMeasureIt, we have developed a multi element SimDis-AED application up to C100+, which measures: Sulfur, Nitrogen and metal distributions. weMeasureIt can be contacted at 510-215-7354 or by e-mail info@wemeasureit.com

Since the transition of the AED, JAS has made 10 substantial improvements to have an enhanced AED system.

To mention a few improvements:

- To ensure a longer lifetime of the AED, a new water drawer with a new water pump is designed.
- JAS has updated the AED software for better performance.
- A new JAS High Temperature transfer line will be introduced soon.

From **February 20th till March 31st** we would like to offer **special Pittcon prices** on the following software upgrades and upgrade packages:

JAS G2350-10120SU	JAS ChemStation AED software upgrade including AED browser from \$ 2,067.00 (list price) for \$ 1,490.00 (Pittcon price)
JAS G2350-10120BU	JAS ChemStation AED bundle including Windows, ChemStation and AED software upgrade, AED browser, Dell PC with 17" Flat Panel Monitor from \$ 4,875.00 (list price) for \$ 3,450.00 (Pittcon price)
JAS 2379ME	JAS-WMI AED Multi Element SimDis software from \$ 8,950.00 (list price) for \$ 6,950.00 (Pittcon price)

JAS will have a booth at the 2004 Pittsburgh Conference in Chicago from March 8 till March 11. Visit us at booth 4287 where you can talk to us about your analytical needs and requirements.

Best Regards,

Joint Analytical Systems (Americas) Inc.
Jan Pijpelink
Sales and Marketing Manager
12 Stow Road, Suite 130
Marlton, New Jersey, 08053
Phone: 856-810-7400, Toll free: 800-343-6015, Fax: 856-810-7404
Email: info@jas-usa.com